# NEOSPLINT

## DESCRIPTION

SplintArt are simple and practical splints created to immobilize limbs. The set includes a complete range of splints for the treatment of sprained or fractured limbs, which can be shaped as needed and once applied according to the needs of the case, offer extraordinary rigidity. Due to their lightness and adaptability, the complete set or the single pieces can be taken everywhere as they take up very little space. Moreover, being made of neoprene, they maintain body temperature.

## COMPOSITION

SplintArt structure in Royal Blue imitation neoprene laminated with MTP mm.5 D90 and black cotton jersey  
Cold mouldable rigid aluminium core  
Blue and orange Velcro closure and adjustment

## USE DESTINATION

Sprained limbs  
Fractured limbs

CODE	MEASURE	Code R.D.M.	Code C.N.D.
STE125	Leg	74683/R	M030599
STE123	Arm	74727/R	M030599
STE124	Elbow - Ankle	74731/R	M030599
STE122	Forearm	74732/R	M030599
STE126	Wrist	74737/R	M030599
STE127	Complete Set (bag + 1 piece to measure)	74761/R	M030599

## METHOD OF APPLICATION

1. Place the injured limb on the SplintArt
  2. Shape the metal core as desired
  3. Adjust and close through Velcro (orange Velcro indicates top)
- (To return the SplintArt to a straight position, straighten the metal core with your hands)

## INDICATIONS

Post-traumatic or post-operative situations of the limbs where there is a need for effective support and/or protection of the skeletal system and muscle bundles.

When using the devices in the presence of skin lesions, protect the lesion with appropriate bandaging; in such situations, do not use the device in direct contact with the lesion.

Sanitary product to be applied under medical supervision.

Radiolucent – Latex Free – Non-Sterile – Single patient.

## CONTRAINDICATIONS / WARNINGS

Disorders of the venous, arterial or lymphatic system. If a medical condition is suspected, consult your physician before use. In any case, use for prolonged periods of time only under medical supervision. It is advisable not to apply the product directly in contact with the skin.

## PRESERVATION

Store in a dry place, away from sunlight, heat sources and dust.

The product has no expiration date, but the expected life span is 10 years for the materials and components of the device.

## CLEANING

-Hand wash in water at 30° C with neutral soap or detergent

-Rinse well

-Do not wring -Let dry horizontally -Do not iron

## DISPOSAL

Not dangerous material, assimilable to urban waste

## PACKAGING

Primary: Pieces 1 in PVC bag

Secondary: box pieces 5

Disposable packaging and packing.

## PRODUCT CONFORMITY/REFERENCE STANDARDS

REGULATION (EU) 2017/745

Class I Medical Device

ISO 9001:2015

PVS SPA  
via Leonardo da Vinci, 18  
20051 Cassina de Pecchi MI IT  
tel 029160011 fax 0291600149 mail info@pvs-spa.it

# NEOSPLINT

## MANUFACTURER

Artsanity S.r.l. – Via della Borsa, 18 – 31033 Castelfranco Veneto (TV) – ITALY

## PACKAGE LEAFLET ACCOMPANYING THE PRODUCT

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#### COMPOSITION

SplintArt structure in Blue Royal neoprene imitation coupled with MTP 5mm D90 and black cotton jersey. Rigid cold moldable aluminum core - Blue and orange velcro closure and adjustment

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#### INSTRUCTION FOR USE

Rest the injured limb on the SplintArt. Model the metal core as desired. Adjust and close through the velcro straps (the orange velcro indicates the upper part) (To return the SplintArt to a straight position, straighten the metal core with the hands)

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#### WARNING

Disorders of the venous, arterial or lymphatic system. In case of suspected ongoing disease, consult your doctor before use. In any case, use for prolonged periods of time only under medical supervision. It is recommended to apply the product not directly in contact with the skin. X-ray transparent - non sterile - single patient - Latex Free  
**“It is advisable to report to Artsanity Srl and the competent national authority any serious accident occurring in relation to the device”**

Product manufactured and distributed in accordance with Regulatory (UE) 2017 / 745

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#### STORAGE AND CLEANING

Store in a cool, dry and clean place away from direct sunlight and heat sources.  
Hand wash in water at 30 °C with soap or neutral detergent. Do not rinse well, do not wring

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