

DATA SHEET

Code/ref.	ACQ414eye wash 500ml ACQ408eye wash 250ml ACQ409eye wash 100ml
Product:	Eye wash
Intended use:	Eye wash and wound irrigation. First aid procedure.
Assembler	PVS spa - Via L. da Vinci 18, 20051 Cassina de Pecchi – Italy
CE marking:	The products defined as medical devices, contained in the first aid kit comply with MDR 745/2017 of 05/04/2017.

Description:

Bottle intended for ocular first aid procedure and wound irrigation.

Dimensions:

(W x D x W): 320 x 410 x 80 mm Weight: 1.20 kg

Kit composition:

2 bottles of NACl 0.9% 500ml - produced by BBraun Melsugen AG - Class IIA Medical Device pursuant to Directive MDD 93/42/EEC

1 polyethylene eye cap/shower, not medical device



^{*} Images are purely indicative



MANAGEMENT AND STORA	GE
The product does not require special storage and conservation con	nditions.
Temperature between +5°C and +25°C	4
Avoid exposing the product to direct sunlight.	*
Keep in a dry place	*
PRECAUTIONS	
There is information to read	\triangle
Read the instructions for use	Ţį.
Do not use if packaging is damaged.	(6)

REGULATORY INFORMATION

The following DROP procedural system, assembled by PVS SpA, complies with the requirements of EU Regulation 2017/745.

PVS has drawn up the Declaration of Conformity for the DROP system in accordance with Article 22 of EU Regulation 2017/74

DISPOSAL AND ELIMINATION INFORMATION

The disposal of materials is indicated on the packaging with specific symbols

ATTACHMENT

- a. BBraun Technical Data Sheet
- b. Declaration of conformity according to ART. 22

Expiry date given by the Manufacturer of the DM