

TECHNICAL SHEET

INFORMATION ABOUT THE MEDICAL DEVICE "DROP EYE WASH STATION"






TRADE NAME	DROP EYE WASH STATION
INTENDED USE	Eye washing and wound irrigation. First aid procedure.
ASSEMBLER	PVS S.p.A – Via L. da Vinci 18, 20051 Cassina de Pecchi – Italy





COMPOSITION OF THE ASSEMBLY

- N. 2 NaCl 0.9% 500ml- produced by BBraun Melsugen AG – Class IIa Medical Device according to MDD Directive 93/42/EEC
- 1 polyethylene eye cap/shower, not medical device
- 1 mirror, no medical device



TECHNICAL SHEET

MANAGEMENT AND STORAGE	
The product does not require special conditions of conservation and storage. Temperature between +5°C and +25°C	
Avoid exposing the product to the direct action of the sun.	
Keep in a dry place	

PRECAUTIONS	
There is information to read	
Read the operating instructions	
Do not use if packaging is damaged	
Expiration date given by the Manufacturer of the DM	

REGULATORY INFORMATION
<p>The following DROP procedural system, assembled by PVS SpA, complies with the requirements of EU Regulation 2017/745.</p> <p>PVS has drawn up the Declaration of Conformity for the DROP system in accordance with Article 22 of EU Regulation 2017/745.</p>

ATTACHMENTS
<ul style="list-style-type: none"> a. DATASHEET BBRAUN b. DECLARATION OF CONFORMITY ACCORDING TO ART. 22

CODE
REF: CS202 DROP EYE WASH STATION



TECHNICAL SHEET

INFORMATION ON DISPOSAL AND DISPOSAL
The disposal of materials is indicated on the package with appropriate symbology